

UBAC2 antibody



Catalog Number: 116472

Product name

UBAC2 antibody

Immunogen

[Human UBAC2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

UBAC2 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HEK-293 cells, mouse kidney tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of UBAC2 recombinant protein (Accession Number: XM_047430286). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116472(UBAC2 Antibody) at dilution of 1:600



Immunofluorescent analysis of HEK-293 cells using Catalog No:116472(UBAC2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)