

Anti-Fyn (4B8) Mouse antibody



Catalog Number: 168472

Product name

Anti-Fyn (4B8) Mouse antibody

Specificity

Human, Monkey, Rat, Mouse

Antibody description

Fyn (4B8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Fyn protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


WB, ICC/IF

Dilutions


WB: 1/500

;ICC: 1/50

Validations



Immunocytochemistry staining of HeLa cells fixed with -20u2103 Methanol and using Fyn mouse mAb (dilution 1:50).



Western blot detection of Fyn in HeLa,3T3,C6,COS7,CEM,Ramos and Jurkat cell lysates using Fyn mouse mAb (1:500 diluted).Predicted band size:59KDa.Observed band size:59KDa.