

COQ7 antibody

Catalog Number: 109472

Product name

COQ7 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

COQ7 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, mouse heart tissue. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of COQ7 recombinant protein (Accession Number: XM_047433494). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109472(COQ7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:109472(COQ7 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded

COQ7 antibody

Catalog Number: 109472

human liver cancer using Catalog No:109472(COQ7 antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-COQ7 (IP:Catalog No:109472, 4ug; Detection:Catalog No:109472 1:500) with mouse liver tissue lysate 6000ug.

