

# GAS6 antibody

Catalog Number: 110872

## Product name

GAS6 antibody

## Specificity

Human; other species not tested.

## Antibody description

GAS6 Rabbit Polyclonal antibody. Positive WB detected in DU 145 cells. Positive IHC detected in human lung tissue, human ovary tissue. Observed molecular weight by Western-blot: 75 kDa

## Preparation

This antibody was obtained by immunization of GAS6 recombinant protein (Accession Number: NM\_000820). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:110872(GAS6 Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human lung using Catalog No:110872(GAS6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung using Catalog No:110872(GAS6 antibody) at dilution of 1:50 (under 40x lens)