

# TBCD antibody

Catalog Number: 115872

## Product name

TBCD antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TBCD Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, MCF-7 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 130 kDa

## Preparation

This antibody was obtained by immunization of TBCD recombinant protein (Accession Number: NM\_005993). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115872(TBCD antibody) at dilution of 1:500



IP Result of anti-TBCD (IP:Catalog No:115872, 4ug; Detection:Catalog No:115872 1:500) with mouse brain tissue lysate 3440ug.