

alpha N Catenin antibody

Catalog Number: 107872

Product name

alpha N Catenin antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

alpha N Catenin Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells. Positive WB detected in SH-SY5Y cells, mouse pancreas tissue. Positive IP detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 102 kDa




Preparation

This antibody was obtained by immunization of alpha N Catenin recombinant protein (Accession Number: NM_001399737). Purification method: Antigen affinity purified.

SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107872(CTNNA2 antibody) at dilution of

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

1:1200

Clonality

Polyclonal

IP Result of anti-CTNNA2 (IP:Catalog No:107872, 4ug; Detection:Catalog No:107872 1:500) with SH-SY5Y cells lysate 1600ug.

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IP

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000



Immunofluorescent analysis of SH-SY5Y cells

alpha N Catenin antibody



Catalog Number: 107872

using Catalog No:107872(CTNNA2 Antibody) at

dilution of 1:50 and Alexa Fluor 488-conjugated
AffiniPure Goat Anti-Rabbit IgG(H+L)