

EIF1 antibody



Catalog Number: 110172

Product name

EIF1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Immunogen

[Human EIF1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EIF1 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of EIF1 recombinant protein (Accession Number: NM_005801). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110172(EIF1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:110172(EIF1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded

EIF1 antibody



Catalog Number: 110172

human prostate cancer tissue slide using Catalog

No:110172(EIF1 Antibody) at dilution of 1:50
(under 40x lens)