

Anti-SCGN / Secretagogin antibody



Catalog Number: 103172

Product name

Anti-SCGN / Secretagogin antibody

Immunogen

[Human SCGN / Secretagogin \(His Tag\) recombinant protein](#)

Specificity

Human SCGN / Secretagogin

Antibody description

Rabbit polyclonal to SCGN / Secretagogin

Preparation

Produced in rabbits immunized with purified, recombinant Human SCGN / Secretagogin (rh SCGN / Secretagogin; NP_008929.2; Asp 2-Pro 276). SCGN / Secretagogin specific IgG was purified by Human SCGN / Secretagogin affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human SCGN. The detection limit for Human SCGN is approximately 0.00975 ng/well.

IHC-P: 0.1-2 µg/mL

Validations



SCGN / Secretagogin Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human SCGN in human pancreas with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).