

# PRDM10 antibody

Catalog Number: 114172

## Product name

PRDM10 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PRDM10 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in HEK-293 cells, mouse testis tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 120-150 kDa

## Preparation

This antibody was obtained by immunization of PRDM10 recombinant protein (Accession Number: NM\_199437). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114172(PRDM10 antibody) at dilution of



1:1000

IP Result of anti-PRDM10 (IP:Catalog No:114172, 5ug; Detection:Catalog No:114172 1:500) with mouse testis tissue lysate 8000ug.