

ZNF395 antibody



Catalog Number: 117172

Product name

ZNF395 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human; other species not tested.

Validations

Antibody description

ZNF395 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in human brain tissue, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 62 kDa, 72 kDa




Preparation

This antibody was obtained by immunization of ZNF395 recombinant protein (Accession Number: NM_018660). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117172(ZNF395 antibody) at dilution of 1:400

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:117172(ZNF395 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:1000



Immunofluorescent analysis of HepG2 cells using

ZNF395 antibody



Catalog Number: 117172

Catalog No:117172(ZNF395 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)