

CDR2 antibody

Catalog Number: 109172

Product name

CDR2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CDR2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, Y79 cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 62 kDa

Preparation

This antibody was obtained by immunization of CDR2 recombinant protein (Accession Number: XM_005276502). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109172(CDR2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109172(CDR2 antibody) at dilution of 1:50 (under 10x lens)