

# Human SERPINE2/PI7/PN1 (His tag) recombinant protein



Catalog Number: 519172

## General Information

### Gene Name Synonym

Glia-derived Nexin; GDN; Peptidase inhibitor 7; PI-7; Protease nexin 1; PN-1; Protease nexin I; Serpin E2; SERPINE2; PI7; PN1

### Protein Construction

A DNA sequence encoding the Human SERPINE2/PI7/PN1 Ser20-Pro397 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human SERPINE2/PI7/PN1 (His tag) recombinant protein