

# FAM156A antibody



Catalog Number: 110572

## Product name

FAM156A antibody

## Specificity

human; other species not tested.

## Antibody description

FAM156A Rabbit Polyclonal antibody. Positive IHC detected in human heart.

## Preparation

This antibody was obtained by immunization of FAM156A recombinant protein (Accession Number: BC000867). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human heart using Catalog No:110572(FAM156A antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:110572(FAM156A antibody) at dilution of 1:50 (under 40x lens)