

Anti-Catalase Rabbit antibody



Catalog Number: 168572

Product name

Anti-Catalase Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P

Antibody description

Catalase Rabbit monoclonal antibody

Dilutions

1:1000

1:20-1:100

Preparation

Antigen: A synthetic peptide of human Catalase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of Catalase in K562, HeLa cell lysates using Catalase Rabbit mAb(1:1000 diluted). Predicted band size: 60kDa. Observed band size: 60kDa.