

GIMAP2 antibody



Catalog Number: 110972

Product name

GIMAP2 antibody

Specificity

human; other species not tested.

Antibody description

GIMAP2 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human kidney tissue. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of GIMAP2 recombinant protein (Accession Number: NM_015660). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110972(GIMAP2 antibody) at dilution of 1:300