

KCNK18 antibody



Catalog Number: 111972

Product name

KCNK18 antibody

Specificity

Human; other species not tested.

Antibody description

KCNK18 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, A375 cells, MCF-7 cells. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of KCNK18 recombinant protein (Accession Number: NM_181840). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:111972(KCNK18 Antibody) at dilution of 1:1000