

HNRNPCL1 Polyclonal Antibody



Catalog Number: 162972

Product name

HNRNPCL1 Polyclonal Antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

IF

Antibody description

Polyclonal antibody to HNRNPCL1


Dilutions

IF 1:50 - 1:200

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human HNRNPCL1 (NP_001013653.1).

Validations



Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Immunofluorescence - HNRNPCL1 Polyclonal Antibody

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Immunofluorescence analysis of U2OS cells using HNRNPCL1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.