

Anti-UCHL3 antibody



Catalog Number: 103972

Product name

Anti-UCHL3 antibody

Immunogen

[Rat UCHL3 \(His Tag\) recombinant protein](#)

Specificity

Rat UCHL3

Antibody description

Rabbit polyclonal to UCHL3

Preparation

Produced in rabbits immunized with purified, recombinant Rat UCHL3 (rR UCHL3; Q91Y78; Glu 2-Ala 230). UCHL3 specific IgG was purified by Rat UCHL3 affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions


WB: 10-20 μ g/mL

ELISA: 0.5-1 μ g/mL

This antibody can be used at 0.5-1 μ g/mL with the appropriate secondary reagents to detect Rat UCHL3. The detection limit for Rat UCHL3 is approximately 0.00975 ng/well.


IP: 1-4 μ g/mg of lysate

Validations



Lanes	A
Sample (whole cell lysate)	Jurkat
Sample Volume (μ g/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended concentration	10-20 μ g/ml
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L) at 1:5000 dilution
Antibody	Developed using Odyssey imaging system
Explanation	Predicted band size: 23 kDa Observed band size: 23 kDa Additional bands at: 55 kDa (We are unsure as to the identity of these extra bands. They are possible non-specific bindings).

UCHL3 Antibody, Rabbit PAb, Antigen Affinity



Lanes	A	B
Sample (whole cell lysate)	K562	Jurkat
Sample quantity	0.5 mg	
IP antibody quantity	2 μ g	
Protein G Agarose	15 μ l of 50% Protein G Agarose	
Gel	13% SDS-PAGE reducing gel	
Primary antibody	ratUCHL3 antibody at 5 μ g/ml	
Secondary antibody	Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution.	

Purified, Western blot

UCHL3 Antibody, Rabbit PAb, Antigen Affinity
Purified, Immunoprecipitation