

Anti-F13B antibody



Catalog Number: 175972

Product name

Anti-F13B antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to F13B

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human f13b. aa range:61-110

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of COLO cells using Factor XIII B Polyclonal Antibody diluted at 1:1000