

Anti-TRAF6 Rabbit antibody



Catalog Number: 166972

Product name

Anti-TRAF6 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TRAF6 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human TRAF6

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/500-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200

Validations



Immunohistochemical analysis of paraffin-embedded human lung carcinoma, using TRAF6



Antibody.

Western blot analysis of TRAF6 expression in Jurkat cell lysate.