

# FMO2-specific antibody



Catalog Number: 106972

## Product name

FMO2-specific antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

FMO2-specific Rabbit Polyclonal antibody.  
Positive WB detected in mouse kidney tissue,  
mouse lung tissue. Observed molecular weight by  
Western-blot: 50kd, 57kd

## Preparation

This antibody was obtained by immunization of  
Peptide (Accession Number: NM\_001460).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse kidney tissue were subjected to SDS PAGE  
followed by western blot with Catalog  
No:110729(FMO2-specific antibody) at dilution of  
1:300