

Elafin/Skalp antibody



Catalog Number: 110272

Product name

Elafin/Skalp antibody

Specificity

Human; other species not tested.

Antibody description

Elafin/Skalp Rabbit Polyclonal antibody. Positive IHC detected in human skin cancer tissue, human oesophagus cancer tissue, human oesophagus tissue. Positive IF detected in HeLa cells.

Preparation

This antibody was obtained by immunization of Elafin/Skalp recombinant protein (Accession Number: NM_002638). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, IHC


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemistry of paraffin-embedded human skin cancer tissue slide using Catalog No:110272(PI3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skin cancer tissue slide using Catalog No:110272(PI3 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:110272(PI3 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)