

MGAT4B Polyclonal Antibody



Catalog Number: 161272

Product name

MGAT4B Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to MGAT4B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 409-548 of human MGAT4B (NP_055090.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - MGAT4B Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using MGAT4B antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.