

PTGR1 antibody

Catalog Number: 114272

Product name

PTGR1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PTGR1 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in human liver tissue, HeLa cells, human colon tissue. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of PTGR1 recombinant protein (Accession Number: NM_012212). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114272(PTGR1 antibody) at dilution of 1:500



IP Result of anti-PTGR1 (IP:Catalog No:114272, 4ug; Detection:Catalog No:114272 1:500) with HepG2 cells lysate 3200ug.