

ZXDA antibody



Catalog Number: 117272

Product name

ZXDA antibody

Specificity

Human; other species not tested.

Antibody description

ZXDA Rabbit Polyclonal antibody. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 85-90 kDa

Preparation

This antibody was obtained by immunization of ZXDA recombinant protein (Accession Number: NM_007156). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:117272(ZXDA Antibody) at dilution of 1:300