

FKBP1B antibody



Catalog Number: 110672

Product name

FKBP1B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FKBP1B Rabbit Polyclonal antibody. Positive WB detected in U-937 cells, Jurkat cells, mouse brain tissue, mouse heart tissue. Observed molecular weight by Western-blot: 12kd

Preparation

This antibody was obtained by immunization of FKBP1B recombinant protein (Accession Number: NM_054033). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:110672(FKBP1B antibody) at dilution of 1:1000