

Anti-CD11a Rabbit antibody



Catalog Number: 167672

Product name

Anti-CD11a Rabbit antibody

Specificity

Mouse

Antibody description

CD11a Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of mouse CD11a

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/500

Validations



Western blot detection of CD11a in A549,HL-60 cell lysates using CD11a Rabbit mAb(1:1000 diluted).Predicted band size:128kDa.Observed band size:160kDa.