

Human CPB1/CPB/PCPB (His tag) recombinant protein



Catalog Number: 519672

General Information

Gene Name Synonym

Carboxypeptidase B, Pancreas-Specific Protein, PASP, CPB1, CPB, PCPB

Protein Construction

A DNA sequence encoding the Human CPB1/CPB/PCPB His16-Tyr417 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM MES, 150mM NaCl, pH 5.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human CPB1/CPB/PCPB (His tag) recombinant protein