

Anti-ATP6V1H antibody



Catalog Number: 176372

Product name

Anti-ATP6V1H antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to ATP6V1H

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human atp6v1h. aa range:341-390

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using V-ATPase H Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000