

Human TNFRSF25/APO3/DDR3 (Fc tag) recombinant protein



Catalog Number: 518372

General Information

Gene Name Synonym

Tumor necrosis factor receptor superfamily member 25,APO3, DDR3, DR3, TNFRSF12, WSL, WSL1,TNFRSF25,Protein WSL-1,LARD,Protein WSL,Lymphocyte-associated receptor of death,Death receptor 3,Apoptosis-mediating receptor TRAMP,Apoptosis-inducing receptor AIR,Apoptosis-mediating receptor DR3,Apo-3

Protein Construction

A DNA sequence encoding the Human TNFRSF25/APO3/DDR3 Gln25-Phe201 is expressed with a Fc tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of PBS,pH7.4.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human TNFRSF25/APO3/DDR3 (Fc tag)
recombinant protein