

# TMCC1 antibody



Catalog Number: 116082

---

## Product name

TMCC1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

TMCC1 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, MCF7 cells, mouse brain tissue. Observed molecular weight by Western-blot: 72kd

## Preparation

This antibody was obtained by immunization of TMCC1 recombinant protein (Accession Number: NM\_001349264). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116082(TMCC1 antibody) at dilution of 1:500