

BCRP/ABCG2 antibody



Catalog Number: 107082

Product name

BCRP/ABCG2 antibody

Specificity

Human; other species not tested.

Antibody description

BCRP/ABCG2 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Positive IP detected in NCCIT cell. Observed molecular weight by Western-blot: 70-80; 140-160 kDa

Preparation

This antibody was obtained by immunization of BCRP/ABCG2 recombinant protein (Accession Number: NM_004827). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107082(BCRP,ABCG2 antibody) at dilution of



IP Result of anti-BCRP,ABCG2 (IP:Catalog No:107082, 4ug; Detection:Catalog No:107082 1:1000) with NCCIT cell lysate 2000ug.