

Anti-Cystatin C (7F6) Mouse antibody



Catalog Number: 168082

Product name

Mouse IgG1

Anti-Cystatin C (7F6) Mouse antibody

Applications

Specificity

sELISA Detector

Human

Validations

Antibody description

Cystatin C (7F6) Mouse monoclonal antibody

Preparation

Antigen: Recombinant human cystatin c protein.


Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Standard Curve for Human Cystatin C: Capture Antibody Mouse mAb [6F12-C7-D8] to Human Cystatin C at 4u03bcg/ml and Detector Antibody Mouse mAb [7F6-A5-F3]to Human Cystatin C at 0.5u03bcg/ml.