

Anti-IFNGR1 Rabbit antibody



Catalog Number: 169082

Product name

Anti-IFNGR1 Rabbit antibody

WB: 1/200-1/1000

IHC: 1/20-1/100

Specificity

Human

Antibody description

IFNGR1 Rabbit polyclonal antibody

Preparation

Antigen: Synthetic peptide of human IFNGR1

Formulation

-20°C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, ELISA

Dilutions

Validations



The image is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 169082(IFNGR1 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original



magnification: ×200)

Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane: HEPG2 cell lysate, Primary antibody: 169082(IFNGR1 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds