

# Anti-CD98/SLC3A2 antibody



Catalog Number: 100882

## Product name

Anti-CD98/SLC3A2 antibody

## Immunogen

[Mouse CD98/SLC3A2 \(His Tag\) recombinant protein](#)

## Specificity

Mouse SLC3A2

**No cross-reactivity** in ELISA with Human SLC3A2

## Antibody description

Rabbit monoclonal to CD98/SLC3A2

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse SLC3A2 (rM SLC3A2; NP\_001154885.1; Ala139-Ala565).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 1-10 µg/mL

## Validations



SLC3A2 / CD98 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse SLC3A2 in mouse kidney with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to nephric tubule.