

GATA4-Specific antibody

Catalog Number: 110882

Product name

GATA4-Specific antibody

IP: 1:200-1:2000

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

GATA4-Specific Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human heart tissue. Positive IP detected in mouse heart tissue. Positive WB detected in mouse liver tissue, mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 45kd




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002052). Purification method: Antigen affinity purified.

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110882(GATA4-Specific antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human testis using Catalog No:110882(GATA4-Specific antibody) at dilution of 1:50 (under 10x lens)




GATA4-Specific antibody

Catalog Number: 110882

Immunohistochemical of paraffin-embedded human testis using Catalog No:110882(GATA4-Specific antibody) at dilution of 1:50 (under 40x

IP Result of anti-GATA4-Specific (IP:Catalog No:110882, 4ug; Detection:Catalog No:110882 1:600) with mouse heart tissue lysate 3200ug.



lens)