

MTUS1 antibody



Catalog Number: 112882

Product name

MTUS1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MTUS1 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, mouse uterus tissue. Observed molecular weight by Western-blot: 120kd,180kd

Preparation

This antibody was obtained by immunization of MTUS1 recombinant protein (Accession Number: BC033842). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112882(MTUS1 antibody) at dilution of 1:500