

# HSF2 Polyclonal Antibody



Catalog Number: 163882

## Product name

HSF2 Polyclonal Antibody

## Specificity

Human, Rat

## Antibody description

Polyclonal antibody to HSF2

## Preparation

Antigen: Recombinant protein of human HSF2

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Western blot - HSF2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HSF2 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - HSF2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using HSF2 antibody at dilution of 1:200 (40x



Immunofluorescence - HSF2 Polyclonal Antibody

# HSF2 Polyclonal Antibody



Catalog Number: 163882

---

Immunofluorescence analysis of A549 cells using HSF2 antibody . Blue: DAPI for nuclear staining.