

ALPPL2 antibody



Catalog Number: 107882

Product name

ALPPL2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ALPPL2 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue, human small intestine tissue.

Preparation

This antibody was obtained by immunization of ALPPL2 recombinant protein (Accession Number: NM_031313). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:107882(ALPPL2 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:107882(ALPPL2 Antibody) at dilution of 1:200 (under 40x lens).