

Anti-SPINK5/LEKTI Rabbit antibody



Catalog Number: 168882

Product name

Anti-SPINK5/LEKTI Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SPINK5/LEKTI Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human SPINK5

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/20-1/100

Validations



Western blot detection of SPINK5/LEKTI in HeLa cell lysates using SPINK5/LEKTI Rabbit mAb(1:500 diluted).Predicted band size:121kDa.Observed band size:121kDa.