

LMOD1 Polyclonal Antibody



Catalog Number: 165182

Product name

LMOD1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to LMOD1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 29-187 of human LMOD1 (NP_036266.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Validations



Immunohistochemistry - LMOD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using LMOD1 antibody at dilution of 1:200 (40x lens).