

# USH1G antibody

Catalog Number: 116582

## Product name

USH1G antibody

## Specificity

Human; other species not tested.

## Antibody description

USH1G Rabbit Polyclonal antibody. Positive IP detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 52-55 kDa

## Preparation

This antibody was obtained by immunization of USH1G recombinant protein (Accession Number: NM\_173477). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:2000

## Validations



IP Result of anti-USH1G (IP:Catalog No:116582, 4ug; Detection:Catalog No:116582 1:500) with NIH/3T3 cells lysate 1200ug.