

BRD7 antibody

Catalog Number: 108582

Product name

BRD7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BRD7 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, L02 cells, mouse cerebellum tissue, mouse lung tissue, rat brain tissue, SH-SY5Y cells. Positive IHC detected in human colon cancer tissue, human cervical cancer tissue. Observed molecular weight by Western-blot: 70-74kd

Preparation

This antibody was obtained by immunization of BRD7 recombinant protein (Accession Number: NM_013263). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117229(BRD7 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:117229(BRD7 Antibody) at dilution of 1:50 (under 40x lens)