

GIT1 antibody



Catalog Number: 110982

Product name

GIT1 antibody

Specificity

Human; other species not tested.

Antibody description

GIT1 Rabbit Polyclonal antibody. Positive WB detected in ROS1728 cells.

Preparation

This antibody was obtained by immunization of GIT1 recombinant protein (Accession Number: BC067358). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



ROS1728 cells were subjected to SDS PAGE followed by western blot with Catalog No:110982(GIT1 Antibody) at dilution of 1:1000