

Rat CD23/FCER2 (Fc Tag) recombinant protein



Catalog Number: 502982

General Information

Gene Name Synonym

Fc receptor, IgE, low affinity II, alpha polypeptide, isoform CRA_a; Protein Fc ϵ 2

Protein Construction

A DNA sequence encoding the rat FCER2 (Q6AZ45) (Glu50-Pro331) was expressed with Fc region of human IgG1 at the N-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu

Molecular Mass

The recombinant rat FCER2 comprises 542 amino acids and predicts a molecular mass of 60.2 kDa. The apparent molecular mass of the recombinant protein is approximately 61 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, PH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat CD23 / FCER2 / FCER2A Protein (Fc Tag)
SDS-PAGE