

Anti-IDO2 antibody



Catalog Number: 105982

Product name

Anti-IDO2 antibody

Immunogen

[Human IDO2 \(His Tag\) recombinant protein](#)

Specificity

Human IDO2

Antibody description

Rabbit polyclonal to IDO2

Preparation

Produced in rabbits immunized with purified, recombinant Human IDO2 (rh IDO2; Q6ZQW0-1; Met1-Gly407). IDO2 specific IgG was purified by Human IDO2 affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 5-10 μ g/ml

ELISA:0.1-0.2 μ g/mL

This antibody can be used at 0.1-0.2 μ g/mL with the appropriate secondary reagents to detect Human IDO2. The detection limit for Human IDO2 is 0.00975 ng/well.

Validations



Items	Lanes	A
Sample (whole cell lysate)		HepG2
Sample Volume (μ g/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		5-10 μ g/ml
Secondary Antibody		Dylight 800-labeled Antibody To Rabbit IgG (H+L) at 1:5000 dilution

IDO2 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot