

XBP1 antibody



Catalog Number: 116982

Product name

XBP1 antibody

Immunogen

[Human XBP1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

XBP1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of XBP1 recombinant protein (Accession Number: NM_005080). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116982(XBP1 Antibody) at dilution of 1:300