

Anti-beta Catenin Rabbit antibody



Catalog Number: 168982

Product name

Rabbit IgG

Anti-beta Catenin Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, IP

Human, Mouse, Rat

Dilutions

Antibody description

1:2000

beta Catenin Rabbit monoclonal antibody

1:50

Preparation

1:50

Antigen: A synthetic peptide of human beta Catenin

1:20

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of beta Catenin in C6,CHO-K1,HeLa cell lysates using beta Catenin Rabbit mAb(1:1000 diluted).Predicted band size:86kDa.Observed band size:86kDa.