

HCK antibody



Catalog Number: 111282

Product name

HCK antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HCK Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer. Positive WB detected in Raji cells, Jurkat cells. Observed molecular weight by Western-blot: 57kd

Preparation

This antibody was obtained by immunization of HCK recombinant protein (Accession Number: NM_002110). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111282(HCK antibody) at dilution of 1:100



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:111282(HCK antibody) at dilution of 1:50 (under 10x lens)