

H3N2 HA (His Tag) recombinant protein



Catalog Number: 503282

General Information

Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Arg 345) of the influenza hemagglutinin (A/Perth/16/2009 (H3N2)) (ACS71642.1), termed as HA1, was fused with a polyhistidine tag at the C-terminus.

Organism

H3N2

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The secreted recombinant HA1 subunit of influenza A H3N2 HA (A/Perth/16/2009(H3N2)) comprises 340 amino acids and has a predicted

molecular mass of 38 kDa. As a result of glycosylation, it migrates as an approximately 60 kDa band in reduced SDS-PAGE.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H3N2 (A/Perth/16/2009)
Hemagglutinin Protein (HA1 Subunit) (His Tag)
SDS-PAGE